MANUALE D'USO

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DESCRIPTION OF THE DEVICES

The **SBI** line is indicated for the treatment of alterations in the morphology of the foot.

INTRODUCTION

Optima Molliter S.r.l. declares that SBI walkers are medical devices conceived and designed following the current state of the art, and have been tested to be in compliance with the Regulation EU 2017/745 in its applicable sections. SBI walkers are in compliance with the International Guidelines of the IWGDF for the management of the diabetic foot ulcerations. <https://iwgdfguidelines.org/guidelines/guidelines>. Nevertheless, bear in mind that inappropriate use of the product or underestimation of resulting dangers can cause injuries or harm. They must be correctly used to be effective. It is important to read this User Manual carefully before using the device and to report any serious incident occurring in relation to the device to the manufacturer and to the competent authority of the State in which the user and/or patient is located.

INTENDED USE

SBI MOTUS Treatment of surgical or other plantar ulcer/lesions. /Diabetic foot ulcer / foot ulcer. / Rehabilitation after minor amputation (except for Lisfranc and Chopart amputations). / Management of Charcot foot during transition from acute to chronic phase. In the acute phase if made irremovable with NRS system as 2nd option of the management of Charcot foot.

SBI FRAME Foot stabilization / Ankle fractures - serious ankle sprains. / Metatarsal fractures. / Postoperative use. Management of Charcot foot during transition from acute to chronic phase. In the acute phase with the option of removability as 3rd option of the management of Charcot foot.

SBI W-HEEL Heel wound/ulcer management. / Management of heel ulcer (plantar and dorsal) and Ulcer of Achilles tendon. / Heel decubitus ulcer (plantar and dorsal). / Postoperative rehabilitation of heel/Achilles tendon.

SBI 3x3 PLANTAR INSOLE Plantar insole for the redistribution of the plantar pressure in case of plantar ulcers of the forefoot, midfoot and heel and for the customized offloading of specific areas, with 27 different combinations of offloading.

SBI PLTM PLANTAR INSOLE Plantar insole used to protect the amputated foot and to improve its position inside the brace.

SIDE EFFECTS AND WARNINGS

No side effects related to the use of the devices were recorded. It is advisable not to use the medical device directly on the skin, always cover the foot with bandages or sock. Frequently check that no hot spots and irritation appear on the skin. Do not use these devices if there is a decreased sensitivity resulting from anaesthetic administration. Never walk without the walker sole. The product should be released and used only on medical prescription or on the advice of a licensed healthcare worker. The 3x3 KIT multilayer insoles (standard on SBI Motus and SBI W-Heel), the monolayer insole (standard on SBI Frame) and the PLTM insole need to be handled exclusively by specialized operators (doctor or licensed healthcare worker). **Always check the integrity of the devices and notify the manufacturer in case of negative feedback.**

- It is highly recommended to wear a contralateral protective shoe on the foot without lesions with the same height of the outsole of the device.

- Pay utmost care when walking if there are liquids or slippery substances on the floor. In the case of subjects with precarious walking, use appropriate medical aids such as crutches, tripods, walkers as indicated by your healthcare professional. Notify the Manufacturer in case of any adverse reaction which has not been mentioned in this manual. In the event of a serious accident, notify the Manufacturer and the Competent Authorities. Keep away from heat sources. To be used by subjects weighing up to 130 Kg.

IMPROPER USE The manufacturer **Optima Molliter S.r.l.** shall not be held liable if the product is used under conditions listed below as improper use. Failure to comply with the therapeutic indications / failure to comply with the rules and instructions of the manufacturer / non-compliance with or use contrary to current safety standards / significant alteration of dimensions and shape / significant alteration of the structure by changing-adding items / use product to practice sports, even if not extensive/ use of unsuitable cleaning products / use the product like a tool or like a toy / unsuitable storage / driving automobiles while wearing SBI products / device components are assembled/disassembled by non-qualified personnel / failure to comply with assembly or disassembly instructions / walking without the SBI outsole.

DURATION OF USE

The product retains its characteristics when properly used and stored over its life cycle (on average one season/ 6 months), according to the environmental characteristics and its use. **SBI 3x3** and **SBI PLTM** insoles are subject to greater wear and therefore their operational life is estimated at around 21 days.

HOW TO APPLY AN SBI DEVICE

Open the straps and remove the front tongue. Put on the protective sock included with the device. Once seated, put your foot into the walker by positioning the heel in the rear section of the SBI device.

Position the tongue (tibial protection **1+6**) over the top of the foot and leg. Pay attention that the word **SBI** is positioned as illustrated in the picture, near the toes.

Close the straps in the order indicated from **1** to **4** and carefully insert them into their slits.

HOW TO APPLY THE NON-REMOVABLE SYSTEM

The indication for the Non-Removability of the SBI Motus Walker must be given only by the specialist prescribing the device after assessing the patient's clinical conditions; the Non-Removable System **NRS** must be applied only by a specialist or a licensed healthcare professional. It is advisable to detach the **AFO** main unit (**2**) from the walker before inserting the **NRS** strap. Each package of the SBI Motus Walker contains 1 single-use **NRS** (**4**); Divide the clamp lengthwise into two parts. Insert the first strap into the corresponding lateral slit of the walker bearing the acronym **NRS** on it.

Insert the second strap into the opposite slit of the walker bearing the acronym **NRS** on it and close the first button at the back.

Close the second one in the front. Use a permanent marker to write the date the device was applied in the specific space on the **NRS** strap. The **NRS** must be removed and replaced in 7 days.

REMOVING THE OUTSOLE

Simultaneously press the two release push-buttons (**3**) with the same hand; then pull the sole downwards with the other hand. Now, both front locking points can be easily released.

FITTING THE OUTSOLE

Place the outsole on a flat surface, insert the connectors of the tip of the **AFO** or **AFO OPEN** main body in the two front locking points (**8**) into the specific slits of the outsole. Now, press **SBI** towards the outsole until two "clicks" can be heard. Make sure that the outsole is firmly fastened to the **SBI** device. Risk of slipping on wet surfaces.

INSTRUCTIONS FOR THE ORTHOPAEDIC TECHNICIAN OR LICENSED HEALTHCARE PROFESSIONAL


Read carefully the instructions before unlock the Screws and the Buckles.

REMOVING AND FITTING SCREWS-IN BUTTON AND BUCKLES OF SBI

Screws (11-12) and Buckles (7) can be released from the main unit (2) and repositioned in the specific seats easily and quickly. However, it is important not to underestimate the risks deriving from incorrectly fitting and removing them. Each walker is supplied with 8 SBI screws-in buttons and 8 nut elements. Remove the entire upper soft part and the sole of SBI. Use a screwdriver (to unscrew or tighten screw 11) and four-tooth keyway (to lock the female internal nut element 12). Turn the screw counterclockwise or clockwise to unscrew or screw.

REMOVING AND REFITTING THE TIBIAL REINFORCEMENT

FROM THE TIBIAL PROTECTION To remove the tibial reinforcement (6) unhook the 3 buttons on the reinforcement. To refit it, just reposition it, paying attention to the direction, ensuring that the section bearing the writing "SBI" is near the toes, then hook the 3 buttons.

CUSTOMISING THE SBI PLANTAR INSOLE KIT 3X3

The modular system KIT 3X3 included in SBI Motus and SBI W-Heel and optional in SBI Frame is designed to redistribute plantar pressures and manage ulcers/injuries of the forefoot, midfoot and heel.

WARNING!!! The 3x3 KIT must be customized only by doctors, nurses, orthopedic technicians or qualified health personnel. **WARNING!!!** The KIT 3X3 must be used only inside SBI Motus (1 kit comes as standard), SBI W-Heel (1 kit comes already as standard), and SBI Frame (optional) devices. Never use KIT 3X3 in regular footwear or other walkers. Optima Molliter S.r.l. shall not be held liable if the KIT 3X3 is used in devices other than the SBI Frame, SBI Motus or SBI W-Heel ones. **WARNING!!!** The KIT 3X3 must be used always together with all its components well connected to each other inside the SBI device.

KIT 3X3 consists of 9 interconnected parts in a Puzzle system forming 3 different small plantar insoles; the three plantar insoles have 3 densities (**RED** = soft, **BEIGE** = medium, **BLUE** = rigid) and they can be combined vertically and horizontally. Plantar pressure can be easily redistributed by combining these insoles. The procedure consists in simple steps; to reduce pressure on an ulcerated/injured area, the KIT 3X3 can be customised by marking and cutting part of the plantar insole.

PATIENTS WEIGHING LESS THAN 90 KG

Mark the ulcerated/injured plantar point on the **BEIGE** insole. **1.** Cover the plantar injury/ulcer with a transparent sterilised film and using a lipstick or non-permanent marker mark the perimeter of the injury/ulcer directly over the transparent film (**A**). **2.** Insert the entire insole KIT 3X3 into the device and make sure that the **BEIGE** insole is in contact with the foot (**BEIGE-BLUE-RED**). **3.** Put the patient's foot into the device, position the tibial tongue and close the 4 straps. Help the patient to walk (10 steps). In this way, the plantar surface will be marked at the injury point. **4.** Now, open the device and remove the entire KIT 3X3. **5.** Using a pen, mark the perimeter of the imprint with a plus margin of about 6 mm around and 8 mm in front (**B**). **6.** Using a scalpel with a 10 mm-blade perforate the perimeter of the mark with the blade at about 45° (**C**). **7.** Assemble the components according to the desired pressure redistribution (see Diagram A) and insert the entire KIT 3X3 into the device. **8.** Make sure that the insole with the hole is in the intermediate segment and the red segment to be in contact with the foot. **9.** To avoid infections, always protect the ulcer/injury with suitable dressing before using the KIT 3X3. **10.** Once the Protective Sock has been put on, the device can be used by the patient.


DIAGRAM A

PATIENTS WEIGHING MORE THAN 90 KG

Mark the ulcerated/injured plantar point on the **BLUE** insole.

1. Cover the plantar injury/ulcer with a transparent sterilised film and using a lipstick or non-permanent marker mark the perimeter of the injury/ulcer directly over the transparent film (**D**). 2. Insert the entire insole KIT 3X3 into the device and make sure that the **BLUE** insole is in contact with the foot (**BLUE-BEIGE-ROSSO**). 3. Put the patient's foot into the device, position the tibial tongue and close the 4 straps. Help the patient to walk (10 steps). In this way, the plantar surface will be marked at the injury point. 4. Now, open the device and remove the entire KIT 3X3. 5. Using a pen, mark the perimeter of the imprint with a plus margin of about 6 mm around and 8 mm in front (**E**). 6. Using a scalpel with a 10 mm-blade perforate on the **BLUE** insole the perimeter of the mark with the blade at about 45° (**F**). 7. Place the **BLUE** drilled part on the corresponding **BEIGE** part. Mark the base of the conical hole with a pen. Drill a conical hole also on the **BEIGE** insole (**G**). 8. Assemble the components according to the desired pressure redistribution (see Diagram B) and insert the entire KIT 3X3 into the device. 9. Make sure that the **BLUE** insole with the hole is in the intermediate segment, the **BEIGE** insole in the segment below and the red segment to be in contact with the foot. 10. To avoid infections, always protect the ulcer/injury with suitable dressing before using the KIT 3X3. 11. Once the Protective Sock has been put on, the device can be used by the patient.




DIAGRAM B

WARNING!!! THE AREAS WITH THE OFF-LOADING HOLE MUST NEVER BE IN CONTACT WITH THE PLANTAR ULCER/INJURY. **WARNING!!!** NEVER USE SBI WITHOUT THE INSOLE INCLUDED AND THE PROTECTIVE SOCK. **WARNING!!!** REPLACE THE ENTIRE KIT 3X3 every 21 days on the prescription of the specialist (in case the ulcer/injury is not healed yet). For further information, request the video with the instructions for use to the following e-mail address: info@optimamolliter.com or scan the QR code

PREPARING THE PLTM PLANTAR INSOLE

1. Remove the adhesive film of the black Velcro band and attach it on the insole internal surface. 2. Attach the orange material with the Velcro section on the Velcro band. 3. Position the foot with minor amputation on the insole being careful that there is at least 1,5 cm between the foot and the orange material. 4. Always protect foot and leg with adequate dressing. Put on the protective sock supplied with each device over the dressing.

WARNING!!! PLTM must be prescribed only by a specialist physician and customised only by physicians, nurses, orthopaedic technicians or qualified healthcare personnel. The PLTM must be used only inside the SBI. Never use the PLTM in regular footwear or other devices. Optima Molliter S.r.l. shall not be held liable if the PLTM is used in devices other than SBI. PLTM is designed for postoperative following toe amputation or transmetatarsal amputation.

CLEANING OF THE PRODUCT

The cleaning of the product must be carried out as follows: Rigid support (2): in the dishwasher (maximum temperature: 45°) or by hand with mild soap and cold-water rinsing. Soft coating (1), Sole (9-10), tibial reinforcement (6): use cold water and clean with mild soap; rinse once the operation is completed. Straps (5), Buckles (7), Sock (13): washing machine, cold water. Screw buttons (11), and Nut elements (12): once removed, they must be cleaned with a soft cloth and dedicated galvanized steel cleaning product (screw buttons) and a dedicated product for galvanized iron (nut elements). Avoid any corrosive liquids such as alcohol, bleach, ammonia, solvents, peroxides or optical bleaching agents. Avoid using abrasive rubbers such as brushes, steel wool, and sand paper, or a knife blade to remove stains. Optima Molliter S.r.l. shall not be held liable for damages caused by abrasive products that attack the materials and by direct exposure to heat sources. Dry by hand with a soft cloth or air dry away from heat sources; avoid direct sunlight. Do not use driers or drying machines/tumble dryers.

DISPOSAL OF THE PRODUCT

The materials used in SBI products are partially environmentally friendly, and therefore, require some specific measures for product disposal. Dispose of product responsibly. Refer to the Italian and local regulations on the disposal of inert products and plastics. All the parts are 100% recyclable.

WARRANTY

The warranty refers only to materials used and to manufacturing defects and protects the shoe in accordance with the law only if recognized as defective by the Company. For the article, which is assumed to be defective, a complaint can only be requested by fulfilling the "Complaints Form" which can be requested to the Company and sent to resi@optimamolliter.com only in the presence of: - intact, unused article kept in its original packaging and joined by proof of purchase which must not be more than 10 days prior to the claim. With the prior authorization of Optima Molliter S.r.l., the item must be sent in the original packaging together with the proof of purchase, exclusively carriage paid and indicating the return authorization number communicated by Optima Molliter on the package. Returns that should arrive without affixing the authorization number clearly visible on the package will be rejected. Items used or shipped without compliance with the above specifications, or valued not covered by warranty, may be returned to the sender only at his charge. In the event of any recognized liability of the Company, the product will be returned by the Company to the reseller or customer.

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Optima Molliter Srl