

TECHNICAL INFORMATION

OPTIMA FREE



REHABILITATION FOOT ORTHOSIS

Treatment of diabetic ulcers/dorsal trauma, post-operative/traumatic phases of foot fractures.

PRODUCT COMPONENTS

- 1- Leg
- 2- Lining
- 3- Single
- 4- Bottom
- 5- Velcro



The Optima line is indicated for the treatment of alterations in the morphology of the foot.

Optima Molliter S.r.l. declares that the Optima tutors are medical devices conceived and designed according to the state of the art, have been checked and made compliant with EU Regulation 2017/745 for the applicable parts. The braces comply with IWGDF's international guidelines for the management of diabetic foot ulcers. <https://iwgdfguidelines.org/guidelines/guidelines>.

THERAPEUTIC INDICATIONS

- Management of dorsal foot ulcers
- Post-operative forefoot
- Post-operative toes
- Post-operative hallux valgus
- Suitable for securing functional bandages.

TECHNICAL FEATURES

Rehabilitation foot orthosis, thermo-mouldable heel counter, breathable and washable upper. Rigid anti-flexion and anti-torsion insole, Fully Rocking sole with movement control and pressure load redistribution.

TECHNICAL CHARACTERISTICS

MATERIALIALS:

UPPER: PU 70% + PES 20% + COT 10%

LINING POLYESTER: 100%

INSOLE INSERT: POLYAMIDE PA6 70% +GLASS FIBER 30%

BOTTOM SOLE: POLYURETHANE

PROTECTIVE SOCK: 90% COTTON, 7% POLYAMIDE, 3% ELASTANE

OPTIMA PRODUCTS ARE LATEX FREE.

CLOSURE: 3 VELCROS

PACKAGING: WEIGHT: GR 205 / MATERIAL: CARTON

OPERATIVE TEMPERATURE AND STORAGE: - 5° | + 45°

PRODUCTION LOT: AS DESCRIBED ON THE LABEL ON THE FRONT OF THE BOX

BASE COMPONENTS: MONOLAYER

AGGIUNTIVI: OPTIMA KIT 3X3 PLANTAR INSOLE /PROTECTIVE SOCK

PRODUCT CLEANING

To clean the product, use warm or cold water and neutral soap; rinse once the operation is completed. Avoid any corrosive liquids such as alcohol, laundry detergents, bleach, ammonia, and solvents. Avoid using abrasive materials such as brushes, steel wool, and sand paper, or a knife blade to remove stains. Optima Molliter S.r.l. shall NOT be held liable for damages caused by abrasive products, by products which attack the materials, by direct exposure to heat sources.




POSSIBILITY OF COMPLETELY CUTTING THE FIRST STRAP

Class I medical device according to Regulation (EU) 2017/745

In accordance with Recommendations 2-5-6-7c of the IWGDF International Guidelines on the Management of the Diabetic Foot. <https://iwgdfguidelines.org/guidelines/guidelines>

TAGLIE

	S 34-36	M 37-39	L 40-42	XL 43-45	XXL 46-48
A MM.	247	267	288	316	336
B MM.	101	108	115	119	131
C MM.	181	183	194	202	208
PESO GR. senza plantare	359	415	461	559	713



Optima Molliter Srl

Company Certified under: ISO 9001:2015 IT346378 - ISO 13485:2016 IT346377

Via Breda 19/21 62012 Civitanova Marche (MC) Italy T +39 0733 801060

info@optimamolliter.com / www.yda.it / www.optimamolliter.com

