

TECHNICAL INFORMATION

PEONIA



TECHNOLOGY

- 1 REMOVABLE FOOTBED
- 2 ADJUSTABLE CLOSURES
- 3 WIDE AND STABILIZING HEEL CENTER
- 4 ABSENCE OF STITCHING IN THE CRITICAL AREAS
- 5 SELF-SHAPING UPPER FOR ALLUX VALGUS AND OTHER DEFORMITIES



TQURA is a collection of functional sneakers designed to offer ideal support to those suffering from mild pain or foot problems. Indicated for well-being and primary prevention, foot protection and cool down. Optima Molliter Srl declares that TQURA footwear are medical devices and can be used for the purpose of containing personalized or custom-made insoles. TQURA is part of the Optima Molliter Health Program. All footwear and devices are tested in the best Institutes and certified. By choosing TQURA, you choose a conscious and wellness-oriented lifestyle. Our vision is to make people free to move in their own world.

THERAPEUTIC INDICATIONS

- Relaxation and well-being.
- Sensitive and/or rheumatic feet.
- Minor deformities (hallux valgus or varus/hammer toes, etc.).
- Prevention and management of nail trauma.
- Post-trauma and/or post-minor surgery rehabilitation of the forefoot.
- Plantar pain (metatarsalgia, heel pain, fasciitis, heel spurs, etc.).
- Foot re-education for biomechanics and posture.
- Polyarthritis, Morton's neuroma, arthrosis, osteoarthritis when used with custom insoles.
- Diabetic foot, no risk, Risk 0, Risk 1* (+ custom insoles), Risk 2* (+ custom insoles); *only for items in the line with a semi-rigid sole.

TECHNICAL CHARACTERISTICS

OUTSOLE: TPR ANTISLIP MATERIAL

MIDSOLE: REMOVABLE INSOLE

INSOLE: SUEDE LEATHER

HIGH HEEL LAST MEASURE: 5 MM

OUTSOLE DROP MEASURE: 10 MM

UPPER: STRETCH NAPPA LEATHER & STRETCH SUEDE LEATHER

LINING: ASSEMBLED FOAM TEXTILE

LACING SYSTEM: ADJUSTABLE VELCRO

AVERAGE SHOE WEIGHT SIZE 37: GR. 200

SIZE - FITTING: 37- FIT 13.5/14 ADJUSTABLE

PRODUCT CLEANING

The cleaning of the product must be carried out using:



- For the upper in fabric, use a cloth slightly dampened with water and neutral soap, rub gently and dry.
- For leathers: special creams (moisturizing-polishing) specific for the type of leather to be treated and compliant and approved by the CEE, coated and subsequently polished with a clean, soft and non-abrasive dry cloth.
- For suede or similar: use only brushes or soft cloths and only with dry leather to remove any stains (mud or similar).
- Make sure to perform a regular hydration of the leather by applying special creams, especially in case of particularly soft leathers such as nappa, deer etc.
- Do not use alcoholic and corrosive liquids, detergents, alcohol, laundry detergents, bleach, ammonia and solvents.
- Do not use abrasive materials and/or products brushes, steel wool, glass paper, blades or similar, to remove stains or residues.
- The footwear MUST NOT be machine washed or washed by being immersed in or put under water.
- To clean the lower lateral surface of the sole use soft sponges dampened in warm water not exceeding 30°C, and neutral soap.
- Do not dry wet or soaked products in front of heat sources. This may irreparably damage the footwear and the warranty will be void. If the product gets accidentally wet, let it dry away from heat sources and in a dry environment.

Optima Molliter Srl shall NOT be held liable for damages caused by abrasive products, by products which attack the materials, by direct exposure to heat sources and by improper use/maintenance.

Class I medical device according to Regulation (EU) 2017/745

SIZES

	35	36	37	38	39	40	41
A	21,76	22,43	23,10	23,77	24,44	25,11	25,78
B	7,40	7,55	7,70	7,85	8,00	8,15	8,30
C	23,80	24,25	24,70	25,15	25,60	26,05	26,50



Optima Molliter Srl

Company Certified under: ISO 9001:2015 IT346378 - ISO 13485:2016 IT346377

Via Breda 19/21 62012 Civitanova Marche (MC) Italy T +39 0733 801060

info@optimamolliter.com / www.yda.it / www.optimamolliter.com

